



Recombinant Marmota monax Interleukin-12 subunit beta (IL12B)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-YP726757MQG |
| Abbreviation | IL12B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q61729 |
| Product Type | Recombinant Protein |
| Immunogen Species | Marmota monax (Woodchuck) |
| Purity | >85% (SDS-PAGE) |
| Sequence | IWELEKNV YVVELDWHPD TPGEKVVLTC DTPEEDGITW TSEQSSEVLG SGKTLTILVK EFEDAGHYTC RRGGEVLSQM LLLHKNEDG IWSTDILKKK EPENKNLVT C EAKNYSGRFT CWWLTAISTD VNFSVKSHRG SSDPQGVTCG EATLSAERVK IEQREYKKYS VQCQEDNACP TAEETLPITV VDAVHKLKY ENYISSFFIR DIIKPDPPKN LKMKPSKTPQ QVEVTWEYPD SWSTPHSYFS LTFSVQVQ GK KKKRSNTLHV DKTSVTVTCQ KGAKVSVQAR DRYYNSSWSE WATMSCP |
| Source | Yeast |
| Target Names | IL12B |
| Protein Names | Recommended name: Interleukin-12 subunit beta Short name= IL-12B Alternative name(s): Cytotoxic lymphocyte maturation factor 40 kDa subunit Short name= CLMF p40 IL-12 subunit p40 |
| Expression Region | 23-327 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |