



Recombinant Mouse Programmed cell death protein 4 (Pdcd4)

Product Code	CSB-MP726762MO
Abbreviation	Pdcd4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61823
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MDIENEQTLN VNPTDPDNLN DSLFSGDEEN AGTTEEIKNEI NGNWISASTI NEARINAKAK RRLRKNSSRD SGRGDSVSDN GSEAVRSGVA VPTSPKGRLL DRRSRSGKGR GLPKKGGAGG KGVWGTPGQV YDVEEVDVKD PNYDDDDQENC VYETVVLPLD ETAFEKTLTP IIQEYFEHGD TNEVAEMLRD LNLGEMKSGV PVLAVSLALE GKASHREMTS KLLSDLCGTV MSTNDVEKSF DKLLKDLPEL ALDTPRAPQL VGQFIARAVG DGILCNTYID SYKGTVDVCVQ ARAALDKATV LLSMSKGGKR KDSVWGS GGG QQPVNHLVKE IDMLLKEYLL SGDISEAEHC LKELEVPHFH HELVYEAIMV VLESTGESAF KMILDLLKSL WKSSTITIDQ MKRGYERIYN EIPDINLDVP HSYSVLERFV EECFQAGIIS KQLRDLCP SR GRKRFVSEGD GGRLKPESY
Source	Mammalian cell
Target Names	Pdcd4
Protein Names	Recommended name: Programmed cell death protein 4 Alternative name(s): Protein MA-3 Topoisomerase-inhibitor suppressed protein
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein localized to the nucleus in proliferating cells. Expression of this gene is modulated by cytokines in natural killer and T cells. The gene product is thought to play a role in apoptosis but the specific role has not yet been determined. Two transcripts encoding different isoforms have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.