



Recombinant Maternal effect lethal protein 47 (mel-47)

Product Code	CSB-EP726786CXX
Abbreviation	mel-47
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q621Q3
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	MASSTAFYL KDAGFHIRNV PKEWNDWNLF HIFQVFGRVG YCRVAGQSND MQLGFVNMLS AADAEVRKN LHDGNLIGDN YSLKVSDHKN IGGALLPMAS ISVQKLLSSP PAKSGPVHLS SSWLPLNKDI EVEVDYLPSS SAAKDVFAL TILRINDAQS NEKYNAMHEK MNSYAQLVAF DSELQIGYDG VFRESRPSIK RVRISATKL YLVDFGKIIN YEKSKCFQIP KVFQTIPTRV SFCGLDGLTW SEQAIPSFND IREVVQKWGA MENATLHAVT CGFAGAINMI NLFCGNSVLA ERLQRKGICE YLPRHQQPRF AYSRDSLLKH TNAGVTAHIS NDADVVKDLL KKIDGVKSML RELDL
Source	E.coli
Target Names	mel-47
Protein Names	Recommended name: Maternal effect lethal protein 47
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.