



Recombinant Burkholderia mallei Aliphatic sulfonates import ATP-binding protein SsuB (ssuB)

Product Code	CSB-MP726864BNS
Abbreviation	ssuB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62K56
Product Type	Recombinant Protein
Immunogen Species	Burkholderia mallei (strain ATCC 23344)
Purity	≥85% (SDS-PAGE)
Sequence	MTGTTLAATY GPISGADLEA ELAQPRIADG DAQDAAVYER DGGARALPFA SGGAPPDGDR ADVRRAGAG DASVRLTRVS KRYGERAVLA DVDLSIGRGS FVSIVGRSGC GKSTLLRLVA ELETSPAGTL VKRGDGGGAL DTRIMYQEAR LLPWKTVLQN VMLGLGRRRAK DDARAVLDEV GLLARANDWP AQLSGGQRQR VALARALVHR PQLLLDEPL GALDALTRIE MHALIERLWR EHRFTALLVT HDVQEAVALA DRVLLIEAGR IAFDQRVPLD RPRARASAAF AALEDRLVQR VLTGSDAAPA APNAAGPEGA SRGRAAPASG LRWAV
Source	Mammalian cell
Target Names	ssuB
Protein Names	Recommended name: Aliphatic sulfonates import ATP-binding protein SsuB EC= 3.6.3.-
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.