



# Recombinant Bacillus cereus GTPase obg (obg)

<b>Product Code</b>	CSB-MP726907BAAD
<b>Abbreviation</b>	obg
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q634A3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bacillus cereus (strain ZK / E33L)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MFVDQVKIYV KGGDGGNGMV AYRREKYVPK GGPAGGDGGK GADVVFIVEE GLRTLMDFRY QRHFKA DRGQ HGMSKGQHGR KSEDLLVKVP PGTVVKDEKT GQILADLVTH GQTAVIAKGG RGGRGNSRFA TATNPAP EIA ENGEPGQERD VILELKV LAD VGLVGFPSVG KSTLLSVVSS ARP KIAEYHF TTIVPNLGVV ETGDNRSFVM ADLPLIEGA HAGVGLGHQF LRHIERTRVI VHVIDMSGLE GRDPYEDYVT INNELKEYNL RLTERPQVVV ANKMDMPDAE ENLQAFKEKV GDEVKIFPIS AVTKQGVRDL LFEVANLIET TPEFPIHEVV DESDTSVMYK FETEGVKFDI TRES DGTFVI SGYDIEKTFK MTD FSRDES V RRFARQMRGM GIDEALRARG AKDGDIVKIL EYEF EFID
<b>Source</b>	Mammalian cell
<b>Target Names</b>	obg
<b>Protein Names</b>	Recommended name: GTPase obg Alternative name(s): GTP-binding protein obg
<b>Expression Region</b>	1-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.