



Recombinant Bacillus cereus Probable endonuclease 4 (nfo)

Product Code	CSB-EP726911BAAD-B
Abbreviation	nfo
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q634Q3
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MLKIGSHVSM SGKKMLLAAS EEA VSYGATT FMIYTGAPQN TRRKPIEELN IEAGRKHMEQ NGIEEIIVHA PYIINVGN TT KPETFQLGVD FLRMEIERTS ALGVAKQIVL HPGAHVGAGA DAGIQQIIKG LNEVLTPDQT VNIALETMAG KGTECGRSFE EIAKIIDGVK YNEKLSVCFD TCHTHDAGYD IVNDFDGV LN EFDKIVGIDR LQVLHINDSK NVRGAGKDRH ENIGFGHIGY KALHHIVHHP QLTHVPKILE TPYVGEDKKD KKPPYKLEIE MLKNGTFDEG LLEKIKAQ
Source	E.coli
Target Names	nfo
Protein Names	Recommended name: Probable endonuclease 4 EC= 3.1.21.2 Alternative name(s): Endodeoxyribonuclease IV Endonuclease IV
Expression Region	1-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.