



# Recombinant *Cricetulus griseus* Dual specificity mitogen-activated protein kinase kinase 1 (MAP2K1)

<b>Product Code</b>	CSB-BP726940DXU
<b>Abbreviation</b>	MAP2K1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q63980
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Cricetulus griseus</i> (Chinese hamster) ( <i>Cricetulus barabensis griseus</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	PKKKPTPIQ LNPTPDGSAV NGTSSAETNL EALQKKLEEL ELEEQQRNRL EAFLTQKQKV GELKDDDFEK ISELGAGNGG VVFKVSHKPS GLVMARKLIH LEIKPAIRNQ IIRELQVLHE CNSPIYVGFY GVFYSDGEIS ICMEHMDGGS LDQVLKKAGR IPEQILGKVS IAVIKGLTYL REKHKIMHRD VKPSNILVNS RGEIKLCDFG VSGQLIDSMA NSFVGTRSYM SPERLQGTHY SVQSDIWSMG LSLVEMAVGR YPIPPDAKE LELLFGCQVE GDAAETPPRP RTPGRPLSSY GMDSRPPMAI FELLDYIVNE PPAKLPSGVF SLEFQDFVNK CLIKNPAERA DLKQLLVHAF IKRSDAEEVD FAGWLCSTIG LNQPSTPTHA ASI
<b>Source</b>	Baculovirus
<b>Target Names</b>	MAP2K1
<b>Protein Names</b>	Recommended name: Dual specificity mitogen-activated protein kinase kinase 1 Short name= MAP kinase kinase 1 Short name= MAPKK 1 EC= 2.7.12.2 Alternative name(s): ERK activator kinase 1 MAPK/ERK kinase 1 Short name=
<b>Expression Region</b>	2-393
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.