



Recombinant Burkholderia pseudomallei Urease accessory protein UreE (ureE)

Product Code	CSB-YP726999BXW
Abbreviation	ureE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63RL2
Product Type	Recombinant Protein
Immunogen Species	Burkholderia pseudomallei (strain K96243)
Purity	≥85% (SDS-PAGE)
Sequence	MRTIDKRIAP NVRLAATLVA RAPALTLAYD ARCKSRLAAT LDTGEDVALV LPRGTVLRDG DVLVADDGAL VRVAAAHEAV LLVRAPDALT LTRAAYHLGN RHTPVEVGAG CLKLEYDPVL ADMLTRLGAT VERASAPFQP EAGAYGGGHR HGHDATFAED YALAQQVFDE HHGHSHSHDH DHDHDHDHDH QHGPCCSHGH HHGHR
Source	Yeast
Target Names	ureE
Protein Names	Recommended name: Urease accessory protein UreE
Expression Region	1-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.