



Recombinant Burkholderia pseudomallei Guanylate kinase (gmk)

Product Code	CSB-BP727001BXW
Abbreviation	gmk
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63RV7
Product Type	Recombinant Protein
Immunogen Species	Burkholderia pseudomallei (strain K96243)
Purity	≥85% (SDS-PAGE)
Sequence	MTDSNRDGAA AHSLHAGVYP GNLFMVVAPS GAGKSTLVNA LLSKDPEIRL SISYTTTRKPR PGEQDGQHYH FTTVEDFRER HARHEFLESA EVHGNYYGTS RVWIEEQMKI GHDVLLLEIDW QGAQQVKKQF RNAVGIFILP PSLAALEERL KKRGQDEPNV ITRRLLAAGS EIAHAAEAQY VVINETFEHA LAELECIVAA TRLRFTSQYA RHAELFVELG IHLPHAE
Source	Baculovirus
Target Names	gmk
Protein Names	Recommended name: Guanylate kinase EC= 2.7.4.8 Alternative name(s): GMP kinase
Expression Region	1-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.