



Recombinant Mouse Forkhead box protein L1 (Foxl1)

Product Code	CSB-EP727367MO-B
Abbreviation	Foxl1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64731
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MSHLFSPPLA ALAASPLLYV YSPERPGLPL AFAPAAALAG PGRVEPPQKP PYSYIALIAM AIQDAPEQRV TLNGIYQFIM DRFPFYHDNR QGWQNSIRHN LSLNECFVKV PREKGRPGKG SYWTLDPRL DMFENGNYRR RKRKPKPAAG SPEAKRTRVE PPESEVGCDV GSPDLATALP TRAPDRSQSP AVGTARPALL PWGPPEPRDP DADLTVQGAG AVASGQLQRP AHHLGSPCLP APSGSPKGSK SKSFSIDSIL AVRPTPASGA EAPGIPKVP GALGSSLLAA SSGLAPPFNA SLVFDAHVQG GFSQLGIPFL SYFPLQVPEA TVLRFH
Source	E.coli
Target Names	Foxl1
Protein Names	Recommended name: Forkhead box protein L1 Alternative name(s): Forkhead-related protein FKHL11 Transcription factor FKH-6
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.