



Recombinant Human Protein FAM154B (FAM154B)

Product Code	CSB-MP727448HU
Abbreviation	FAM154B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q658L1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGAKSMRSWC LCQICSCGSD YCPYEIVKQP RHVPEEYKPK QGKIDLGTTY KRDLSYKVVQ PVAIVRPLER QVKKGKLDTV PTYKDDYRAW DLHKSELYKP EQTYHPPTVK FGNSTTFQDD FVPQEIKPRQ SFKPSSVVKR STAPFNGITS HRLDYIPHQL ELKFERPKEV YKPTDQRFED LTTHRCDFQG LIGETAKLCR PVHTRVTQNA LFEGSTEFRE SFQPWEIPPP EVKKVPEYVP PTGSMLLNST SHLDYVPYQA NHVVPIRPVS QKRSNNFPFQ GKSIMKEDFP AWESCRQGLI KKQQQIPNPS GKFDGLSTFR SHYVPHELIP TESCKPLNIA FKSSVPFDDV TMYSVEYTPK RQEICPASYP SPPGYIFDNT NSQGHKFFRK IIPAVKAF
Source	Mammalian cell
Target Names	SAXO2
Protein Names	Recommended name: Protein FAM154B
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.