



Recombinant Bacillus licheniformis Protein nrdI (nrdI)

Product Code	CSB-MP727488BQU
Abbreviation	nrdI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65JA6
Product Type	Recombinant Protein
Immunogen Species	Bacillus licheniformis (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
Purity	≥85% (SDS-PAGE)
Sequence	MIQIVFDSKT GNVQRFVDKT PFRNKRKVST EEYLDEPFVL ITYTTGFGEV PKTTEMFLEK NAHLLLGVAA SGNRVWGDNF AKSAEKISKQ YQVPILGKFE LSGTAKDVEL FTQEVERVVT KSSAKMDPVK Q
Source	Mammalian cell
Target Names	nrdI
Protein Names	Recommended name: Protein nrdI
Expression Region	1-131
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.