



# Recombinant *Bacillus licheniformis* Ribosomal RNA small subunit methyltransferase A (rsmA)

<b>Product Code</b>	CSB-MP727520BQU
<b>Abbreviation</b>	rsmA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q65PH9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Bacillus licheniformis</i> (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNKDIATPIR TKEILNKYGF SFKKSLGQNF LIDTNILDRI VDHAGVTER GVIEIGPGIG ALTEQLAKRA KKVTAFEIDQ RLLPILEDTL SPYDNVTVIH QDVLKADVRA VMDEQFQDCD EVMVVANLPY YVTTPIIMKL LEENLPLKGI VVMLQKEVAD RMAAKPSSKE YGSLSIQVQF YTEAKTVMNV PKTVFVPQPN VDSAVIRLTL RKEPAVAVQD AAFFQVVK A SFAQRRKTLF NNLVNNLPNG KENKSKIERA LQDSHIDGKR RGESLSIEEF AVLSDRLEVV LL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rsmA
<b>Protein Names</b>	Recommended name: Ribosomal RNA small subunit methyltransferase A EC=2.1.1.182 Alternative name(s): 16S rRNA (adenine(1518)-N(6)/adenine(1519)-N(6))-dimethyltransferase 16S rRNA dimethyladenosine transferase 16S rRNA dimethylase
<b>Expression Region</b>	1-292
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.