



Recombinant Mannheimia succiniciproducens

Amino-acid acetyltransferase (argA)

Product Code	CSB-MP727561MAAB
Abbreviation	argA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65UJ9
Product Type	Recombinant Protein
Immunogen Species	Mannheimia succiniciproducens (strain MBEL55E)
Purity	≥85% (SDS-PAGE)
Sequence	MRSTELVQWF RQSTPYVNMH RGKTFVIMLD GNTIASSNFI NIINDISLLH SLGIKLIIVY GARVQINSLL AQNNVTSVYH KNIRVTDPR T LELVKQAVGQ LSYDITARLS VRLPHSPVLN VVSSNFILAQ PIGVDDGVDY MLSGKIRRIE IDNIKHLDN NAIVLLGPIA PSVTGETFNL PFEEIATQVA IKLKAELIG FSSTQGILDP QGISIPDLLP QDAAKYLNQY IQQGEYHCSQ ARFLQAAIEV CKAGVKRSHL LSYEEDGSL QELFTRDGVG TQLSVDNSED IRIATVQDIP GLIELIHPLE QQGILVKRSR EQLEMDIAN Y TIIDRDGVII ACAALNQYPE ENMAEMACVA VHPDYRSSSR GDILLEAIQK RARQLGIEKL FVLTTTRTVHW FQERGFR LAN VEDLPKEKRD HYN YQRRSKI LIQPLNEEE
Source	Mammalian cell
Target Names	argA
Protein Names	Recommended name: Amino-acid acetyltransferase EC= 2.3.1.1 Alternative name(s): N-acetylglutamate synthase Short name= AGS Short name= NAGS
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.