



# Recombinant Echovirus 16 Genome polyprotein

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP727640ECA                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |
| <b>Uniprot No.</b>       | Q66790                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Echovirus 16 (strain Harrington)                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | GAQVSTQKT GAHETLLEAA QGATINYTNI NYYKDAASNS ANRQDFSQDP<br>SKFTEPVKDI MIKSMPALNS PSAEECGYSD RVRISITLGNs TITTQESANV<br>VVAYGRWPKY LEDDQATAED QPTQPDVATC RFYTLESVQW<br>EANSAGWWWK FPEALKDMGL FGQNMYYHYL GRAGYTIHVQ<br>CNASKFHQGC LLVVCVPEAE MGCAKPDENV DATNLTNGEN TCELTAGAAP<br>AEKGKVQTAV CNATMGVAVG NLTIFPHQWI NLRTNNCATI VMPYINSVPM<br>DNMFRHYNFT LMVIPFVPLT SMGGSTYVPI TVTIAPMCAE YNGLRLSTQH Q              |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Protein Names</b>     | Recommended name: Genome polyprotein Cleaved into the following 6 chains:<br>1. Protein VP0 Alternative name(s): VP4-VP2 Protein VP4 Alternative name(s): P1A Virion protein 4 Protein VP2 Alternative name(s): P1B Virion pro                                                                                                                                                                              |
| <b>Expression Region</b> | 2-331                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |