



# Recombinant Serine hydroxymethyltransferase (glyA)

<b>Product Code</b>	CSB-BP727648YAH
<b>Abbreviation</b>	glyA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q667X1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Yersinia pseudotuberculosis serotype I (strain IP32953)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLKREMNIAD YDADLWRAME QEVVRQEEHI ELIASENYTS PRVMQAQGSQ LTNKYAEGYP GKRYGGCEY VDVVEQLAID RAKALFGADY ANVQPHSGSQ ANVAVYSALL KPGDTVLGMN LAHGGHLTHG SPVNFSGKLY NIVPYGIDES GQIDYEDLAR QAEIHKPKMI IGGFSAYSGI VDWAKMREIA DSIDAWFFVD MAHVAGLVAA GVYPNPVPHA HIVTTTTHTKT LAGPRGGLIL AKGGDEDLYK KLNSSVFPGN QGGPLMHVIA GKAAVALKEAM EPEFKIYQQQ VAKNAKAMVA VFLERGYKVV SGGTDNHLFL LDLVDKDITG KDADAALGRA NITVNKNSVP NDPKSPFVTS GVRIGSPAIT RRGFKEAESR ELAGWMCDVL DNINDEATIE RVKQKVLAIC ARLPVYA
<b>Source</b>	Baculovirus
<b>Target Names</b>	glyA
<b>Protein Names</b>	Recommended name: Serine hydroxymethyltransferase Short name= SHMT Short name= Serine methylase EC= 2.1.2.1
<b>Expression Region</b>	1-417
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.