



Recombinant Betaine aldehyde dehydrogenase (betB)

Product Code	CSB-MP727682YAH
Abbreviation	betB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66D53
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MSRYGLQKLY INGAYTDSTS GDTFDAVNPA NGECIAQLQA ANAQDVDKAV AAAKQGQPVW AAMTAMERSR ILRRAVDILR DRNDELAAIE TADTGKPLSE TRSVDIVTGA DVLEYYAGLI PALEGQQIPL RGSFAVYTRR EPLGVVAGIG AWNYPLQIAL WKSAPALAAG NAMIFKPSEV TSLTALKLAG IYTEAGLPAG VFNVLTGSGD QVGQMLTEHP GIAKVSFTGG IASGKKVMAN AAGSTLKDVT MELGGKSPLI IFADADLDKA ADIAMMANFY SSGQVCTNGT RVFVPLALQA AFEQKIVERV KRIHIGDPSD ERTNFGPLVS FQHRDSVMRY IDSGKREGAT LLIGGYSLTE GALAHGAYVA PTVFTHCRDD MQIVREEIFG PVMSILSYQS EEEVIRRAND TEYGLAAGVV TQDLNRAHRV IHQLQAGICW INTWGESAPE MPVGGYKHSG VGRENGISTL EHYTQIKSIQ VELSSFNSVF
Source	Mammalian cell
Target Names	betB
Protein Names	Recommended name: Betaine aldehyde dehydrogenase Short name= BADH EC= 1.2.1.8
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.