



Recombinant Phosphoadenosine phosphosulfate reductase (cysH)

Product Code	CSB-EP727694YAH
Abbreviation	cysH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66ED2
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MSQFNLSELN ALPKAKQAAA LVLVNGQLEH LTAQERVSWA LDNLPGEFVL SSSFGIQA AV CLHLVTRQRP DIPVILTDG YLFPETYRFI DDLTEKLQLN LQVFRAAHSP AWQEARYGKL WEQGVEGIER YNTLNKVEPM NRALEALGAQ TWFAGLRREQ SGGRSQLPVL ALQRGIFKLL PIIDWDNRQV YQYLTQHGLS YHPLWEQGYL SVGDTHHTRK WEPGMSEET RFFGLKRECG LHEG
Source	E.coli
Target Names	cysH
Protein Names	Recommended name: Phosphoadenosine phosphosulfate reductase EC=1.8.4.8 Alternative name(s): 3'-phosphoadenylylsulfate reductase PAPS reductase, thioredoxin dependent PAPS sulfotransferase PAdoPS reductase
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.