



Recombinant Homoserine kinase (thrB)

Product Code	CSB-MP727702YAH
Abbreviation	thrB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66ET8
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MVKIYAPASI GNVSVGFDVL GAAVSPIDGT LLGDCVSVTA AERFSLHNEG RFVSKLPDDP KQNIVYQCWE RFCQEMGKEI PVAMVLEKNM PIGSGLGSSA CSVVAGLMAM NEFCGQPLDK VTLLGMMGEL EGRVSGSIHF DNVAPCYLGG MQLILEQEGY ISQDVPGFSD WLWVMAYPGI KVSTAEARAI LPAQYRRQDC ITHGRNLAGF IHACHTQQPD LAAKMMKDVI AEPYRTQLLP GFAAARQAAQ DIGALACGIS GSGPTLFAVC NDQATAQRMA GWLQNHYLQN DEGFVHICRL DTAGARLLG
Source	Mammalian cell
Target Names	thrB
Protein Names	Recommended name: Homoserine kinase Short name= HK Short name= HSK EC= 2.7.1.39
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.