



Recombinant Danio rerio Coiled-coil domain-containing protein 130 homolog (ccdc130)

Product Code	CSB-MP727732DIL
Abbreviation	ccdc130
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66I85
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MGERKGVNKW YPPDFDPAKH GSINGYYKTH PLRERARKLS QGILIIRFEM PYNIWCDGCK NHIGMGVRYN AEKKKVGNY TPIYRFRMK CHLCVNYIEM QTDPATCDYV IVSGAQRKEE RWDMAENEQI LITTERNEKEK LETDAMYKLD HGGKDKEKLR AAIPSLNELQ EHQSGWKDDF QLNSALRRKF RTEKKVIAEE EEKDNAVRLR TGLSIPLVPE REEDKKLASL LTFQSPDSYE DKKQWKRQEI SSRSWFNSPS SAAGGAAGSL LQKLGQQGRG AAVAKALSSS TSTLPILVRR KSESSKSETN NTMSTILPVD THPAAKDTDA TDLPSIVNNI NVSINTDINS CTTDACKASS SSEEENSIDS CATGKSLVAD YSDSDSGSEV
Source	Mammalian cell
Target Names	ccdc130
Protein Names	Recommended name: Coiled-coil domain-containing protein 130 homolog
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.