



Recombinant *Xenopus laevis* Dynactin subunit 2-A (dctn2-a)

Product Code	CSB-EP727736XBE-B
Abbreviation	dctn2-a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66J30
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MADPKYADLP GIARNEPDVY ETSDLPEDDQ AEFDAAEELTS TSVEHIIVNP NAAYDKFKDK KVGTRSLDFS DRITKSKRTG YESGEYEILG EGIGMKETPQ QKYQRLLHEV QELTQEVEKT QSTVKESAAE EKLTPVALAK QVASLKQQLV STHLEKLLGP DAAINLTDPD GALAKRLLTQ LDAAKTRKNP EGKSPAKGPG PDTENLVTYE LHCRPEQNKF SQAAKMAELE KRLGELEVAV RCDQDTQNPL TVGLQGSCLM DTVEILQAKV NLLDVASLDQ VEARLQSVLG KMNEIAKHKA AIEDADTESK VHQLYETVQK WDSMSGTLPQ VVQRLLTLKQ LHEQAMQFGQ LLTHLDTTQQ MIANSKLDNT NALAMVQKAM KENLATVEDN FSNIDGRIKK LSK
Source	E.coli
Target Names	dctn2-a
Protein Names	Recommended name: Dynactin subunit 2-A
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.