



Recombinant Influenza A virus Matrix protein 1 (M)

Product Code	CSB-BP727761IGT
Abbreviation	M
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67185
Product Type	Recombinant Protein
Immunogen Species	Influenza A virus (strain A/Swine/Hong Kong/127/1982 H3N2)
Purity	≥85% (SDS-PAGE)
Sequence	MSLLTEVETY VLSIVPSGPL KAEIAQRLED VFAGKNTDLE ALMEWLKTRP ILSPLTKGIL GFVFTLTVPS ERGLQRRRFV QNALNGNGDP NNMDRAVKLY RKLKREITFH GAKEVALSYS TGALASCMGL IYNRMGTVTT EVAFGLVCAT CEQIADSQHR SHRQMVTTTN PLIRHENRMV LASHTAKAME QMAGSSEQAA ERMEVASQAR QMVQAMRTIG THPSSSAGLK DDLENLQAY QKRMGVQMQR FK
Source	Baculovirus
Target Names	M
Protein Names	Recommended name: Matrix protein 1 Short name= M1
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.