



Recombinant Rice hoja blanca virus Suppressor of RNA silencing p3 (p3)

Product Code	CSB-EP727775RBAX
Abbreviation	p3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67897
Product Type	Recombinant Protein
Immunogen Species	Rice hoja blanca virus (strain cr) (RHBV) (Rice hoja blanca virus (strain Costa Rica))
Purity	>85% (SDS-PAGE)
Sequence	MNVSLYYSGT PVSSHLLSK NGLSNIVLTC KDLPIPIDLL SLFFDILNER HPSFDEHMFL QMIRKPDDPE NLSVFLKSAI WMLSHKRDLP GHYRLPLTCL VSTYSEYFVE LKPRQPSTKC WFCKIAKDGL PFRVEGVHGF PSEAELYIVP SKEHAIESFE VLSGKKLYRS PSKHHGYLI ASNKPLTSK YVEYDPSKPD TKP
Source	E.coli
Target Names	p3
Protein Names	Recommended name: Suppressor of RNA silencing p3 Alternative name(s): Protein NS3 Protein p3
Expression Region	1-203
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.