



Recombinant Elongation factor Tu (tuf)

Product Code	CSB-MP727797SRW
Abbreviation	tuf
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67JU1
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MAKQRFERTK PHVNIGTIGH VDHGKTTTLTA AITKVLAEKG KAQFMAYDAI DKAPEERERG ITINTSHVEY ETDARHYAHV DCPGHADYVK NMITGAAQMD GAILVVSAAAD GPMPQTREHI LLARQVGVPA IVVFLNKCDM VDDEELLELV ELEVRELLNQ YEFPGDEIPF IRGSALKALE GDPKYVKAIE ELMDAVIDSYI PTPQRDADKP FLMPIEDVFT ITGRGTVVVG RVERGKCKVG DQVEIVGLRE ESKTTVVVTGL EMFRKILDEV QAGDNVIGALL RGIEKKEVER GQVLAKPGSI KPHTKFSGAI YVLTKEEGGR HSPFFNGYRP QFYFRTTDDVT GTIKLPEGVE MVMPGDNVEI AVELIHPIAI EEGLRFAVRE GGRTVASGRV TKVSE
Source	Mammalian cell
Target Names	tuf
Protein Names	Recommended name: Elongation factor Tu Short name= EF-Tu
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.