



# Recombinant Argininosuccinate synthase (argG)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP727802SRW                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Abbreviation</b>      | argG                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                             |
| <b>Uniprot No.</b>       | Q67KE1                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Immunogen Species</b> | Symbiobacterium thermophilum (strain T / IAM 14863)                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Sequence</b>          | MAKKKCVLAY SGGLDTSVAI HWIKENYDAD VVALAVDVGA NDKDLEFIRQ<br>KALSIGAVKS YVYDARKEFA YDFILPTLQA NACYESKYYL SAALSRPLIA<br>KLLVQIAEEE GADFVAHGCT GKGNDQVRFE VSVAALNPNL QVLAPVREWG<br>WSREEEIDYA QKHGIPVPVG KENPFSIDVN LWGRSCEAGV LEDPWAEAPE<br>EAFEWTVNPK DAPDEPEYVE IGFEQGVVPS LNGEALDLVT LIERLHAIAG<br>RHGVGRIDHV ENRLVGIKSR EVYEMPAAAV LTLAHKELET ITLPRELHHF<br>KLGLEQKYAE LVYNGLWFHP LKEALDAFIK KSQETVTGTV RVRLFKGSAF<br>CVGRTSPYTL YSYDLATYDK ADKFDHKAAC GFIDIFGLPT KVYAAVQKSK A |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Target Names</b>      | argG                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Protein Names</b>     | Recommended name: Argininosuccinate synthase EC= 6.3.4.5 Alternative name(s): Citrulline--aspartate ligase                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Expression Region</b> | 1-401                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                    |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                             |