



Recombinant Mouse Transcription factor HES-1 (Hes1)

| | |
|--------------------------|---|
| Product Code | CSB-EP010307MO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P35428 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MPADIMEKNS SSPVAATPAS VNTTPDKPKT ASEHRKSSKP IMEKRRRARI NESLSQLKTL ILDALKKDSS RHSKLEKADI LEMTVKHLRN LQRAQMTAAL STDPSVLGKY RAGFSECMNE VTRFLSTCEG VNTEVRTRLL GHLANCMTQI NAMTYPGQAH PALQAPPPPP PSGPAGPQHA PFAPPPPLV PIPGGAAPPP GSAPCKLGSQ AGEAAKVFGG FQVVPAPDGQ FAFLIPNGAF AHSGPVIPVY TSNSGTSVGP NAVSPSSGSS LTSDSMWRPW RN |
| Source | E.coli |
| Target Names | Hes1 |
| Protein Names | Recommended name: Transcription factor HES-1 Alternative name(s): Hairy and enhancer of split 1 |
| Expression Region | 1-282 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Target Details | This protein belongs to the basic helix-loop-helix family of transcription factors. It is a transcriptional repressor of genes that require a bHLH protein for their transcription. The protein has a particular type of basic domain that contains a helix interrupting protein that binds to the N-box rather than the canonical E-box. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |