



Recombinant Arabidopsis thaliana Zinc finger AN1 domain-containing stress-associated protein 12 (SAP12)

Product Code	CSB-YP727842DOA
Abbreviation	SAP12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67YE6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MAGGGTEAFP DLGEHCQDPD CKLLDFLPFT CDGCKLVFCL EHRYSKSHNC PKSDHGSRTV SICETCSIAI ETTGFDEKGI KSLLEKHERS GDCDPNKKKK PTCPVKRCKE ILTFANNLTC KYCGVKFCLK HRFPTDHVCN KKIINTAGTS SRWNERFMEA LSLRNQKGCG RGSSVSSKSS PSVRSF
Source	Yeast
Target Names	SAP12
Protein Names	Recommended name: Zinc finger AN1 domain-containing stress-associated protein 12 Short name= AtSAP12
Expression Region	1-186
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.