



Recombinant Mouse Putative adenosylhomocysteinase 3 (Ahcyl2)

Product Code	CSB-BP727869MO
Abbreviation	Ahcyl2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68FL4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SVQVVSAAA AAKVPEVELK DLSPSEAEPQ LGLSAAVGA MVPPAGGGDP EAPAPAPAAE RPPAPGPGSG PTAALSPAAG KVPQASAMKR SDPHHQHQRH RDGGEALVSP DGTVTEAPRT VKKQIQFADQ KQEFNKRPTK IGRRSLRSI SQSSTDSSYSS AASYTDSSDD ETSPRDKQQK NSKGSSDFCV KNIKQAEFGR REIEIAEQEM PALMALRKRA QGEKPLAGAK IVGCTHITAQ TAVLMETLGA LGAQCRWAAC NIYSTLNEVA AALAESGFPV FAWKGESEDD FWWCIDRCVN VEGWQPNMIL DDGGDLTHWI YKKYPNMFKK IKGIVEESVT GVHRLYQLSK AGKLCVPAMN VNDSVTKQKF DNLYCCRESI LDGLKRTTDM MFGGKQVVVC GYGEVGKGCC AALKAMGSIV YVTEIDPICA LQACMDGFRL VKLNEVIRQV DIVITCTGNK NVVTREHLDR MKNSCIVCNM GHSNTEIDVA SLRTPELTWE RVRSQVDHVI WPDGKRIVLL AEGRLNLSC STVPTFVLSI TATTQALALI ELYNAPEGRY KQDVYLLPKK MDEYVASLHL PTFDAHLTEL TDEQAKYLGL NKNGPFPKPNY YRY
Source	Baculovirus
Target Names	Ahcyl2
Protein Names	Recommended name: Putative adenosylhomocysteinase 3 Short name= AdoHcyase 3 EC= 3.3.1.1 Alternative name(s): S-adenosyl-L-homocysteine hydrolase 3 S-adenosylhomocysteine hydrolase-like protein 2
Expression Region	2-613
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.