



Recombinant Human Dual specificity phosphatase DUPD1 (DUPD1)

Product Code	CSB-MP727876HU
Abbreviation	DUPD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68J44
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTSGEVKTSL KNAYSSAKRL SPKMEEEGEE EDYCTPGAFE LERLFWKGGSP QYTHVNEVWP KLYIGDEATA LDRYRLQKAG FTHVLNAAHG RWNVDTGPDY YRDMDIQYHG VEADDLPTFD LSVFFYPAAA FIDRALSDDH SKILVHCVMG RSRSATLVLA YLMIHKDRTL VDAIQQVAKN RCVLPNRRGFL KQLRELDKQL VQRRRRSQRQ DGEEEDGREL
Source	Mammalian cell
Target Names	DUPD1
Protein Names	Recommended name: Dual specificity phosphatase DUPD1 EC= 3.1.3.16 EC= 3.1.3.48 Alternative name(s): Dual specificity phosphatase 27
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.