



Recombinant Rickettsia typhi Ribosomal RNA large subunit methyltransferase E (rlmE)

Product Code	CSB-MP727926RNE
Abbreviation	rlmE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68XK5
Product Type	Recombinant Protein
Immunogen Species	Rickettsia typhi (strain ATCC VR-144 / Wilmington)
Purity	≥85% (SDS-PAGE)
Sequence	MTNNLSGYRN KIVRVKTAKK RPISSSNWLR RQLNDPYVAK ARLEGFRSRA AYKLEIHDK FRLFTP NMKI VDLGAAPGGW SQVASKLIKA SDNNLNKII SIDLLKIEPI VGVEFLQKDF FEKDTEELII HGLDGKADLV MSDMASNTIG HKATDHIRTLLLCEQAFEFA LKVLKPSCHF IAKIFRGGAE TALLNKVKRE FRTVKHF KPVS SRSESTEIY LVALNKK
Source	Mammalian cell
Target Names	rlmE
Protein Names	Recommended name: Ribosomal RNA large subunit methyltransferase E EC=2.1.1.166 Alternative name(s): 23S rRNA Um2552 methyltransferase rRNA (uridine-2'-O-)-methyltransferase
Expression Region	1-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.