



Recombinant Gloydus brevicaudus Disintegrin adinbitor

Product Code	CSB-MP727949GGO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q698K8
Product Type	Recombinant Protein
Immunogen Species	Gloydus brevicaudus (Korean slamosa snake) (Agkistrodon halys brevicaudus)
Purity	≥85% (SDS-PAGE)
Sequence	TP AHQRYIELVI VADHGMFTKY NGDSKIREW VRQMVNTVDE IYSYMYIDVA LAGLEIWSNE DLINVQPAAP HTLDSFGKWR ERDLLHRIHH DNAMLLTAID FDGPTIGLAY VGTMCCKPKGS TGVVQDHSTI NLRVAVTMAH EIGHNLGIHH DTGSCSCGGY SCIMSPVISH EPSKYFSDCS YTQCWDFIMN QKPQCILNKP
Source	Mammalian cell
Protein Names	Recommended name: Disintegrin adinbitor
Expression Region	29-230
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.