



Recombinant *Oryza sativa* subsp. *japonica* GTP-binding nuclear protein Ran-3 (RAN3)

Product Code	CSB-BP727964OFG
Abbreviation	RAN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69XM7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSRAQALPDP AAVGYPSFKL ILVGDGGTGK TTFVKRHITG EFEKRYEPTI GVEVRPLDFH TSRGKVRFCWDTAGQEKFGLRDLGYYIHG HCAIIMFDVT SRLTYKNVPT WHKDICRVCD NIPIVLCGNK VDMKNRQVKA KMVTFHRKKN LQYYEISAKS NYNFEKPFLY LARKLTGDMN LRFVEELALL PADVTIDLIA QQKIETEIAA AAAMPLPDED EDGLMD
Source	Baculovirus
Target Names	RAN3
Protein Names	Recommended name: GTP-binding nuclear protein Ran-3 Short name= OsRan3 Alternative name(s): Ras-related nuclear protein 3
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.