



# Recombinant Mouse Kazrin (Kazn)

<b>Product Code</b>	CSB-MP727970MO
<b>Abbreviation</b>	Kazn
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q69ZS8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MMEDNKQLAL RIDGAVQSAS QEVTNLRAEL TATNRRLAEL SGGGGGPGSG PGAATSASAA AVTVADSAVA TMENHQHGAQ VLLREEVVQL QEEVHLLRQM KEMLAKDLEE SQGGKCSEVL SATELRVQLV QKEQELARAK EALQAMKADR KRLKGEKTDL VSQMQQLYAT LESREEQLRD FIRNYEQHRK ESEDAVKALA KEKDLLEREK WELRRQAKEA TDHAAALRSQ LDLKDNRMKE LEAELAMAKQ SLATLTKDVP KRHSLAMPGE TVLNGNQEWV VQADLPLTAA IRQSQQTLYH SHPPHPADRQ VRVSPCHSRQ PSVISDASAA EGDRSSTPSD INSPRHRTHS LCNGDSPGPV QKSLHNPIVQ SLEDLEDQKR KKKKEKMGFG SISRVFARGK QRKSLDPGLF DDSDSQCSPT RHSLSLSEGE EQMDRLQHVE LVRTTPMSHW KAGTVQAWLE VVMAMPMYVK ACAENVKSGK VLLSLSEDL ELGLGVCSSL HRRKLRLAIE DYRDAEAGRS LSKAADLDHH WVAKAWLNDI GLSQYSQAFQ NHLVDGRMLN SLMKRDLEKH LNVSKKFHQV SILLGIELLY QVNFREALQ ERRARCETQN TDPVVWTNQR VLKVVWRDIDL KEYADNLNNS GVBGAVLVLE PTFNAEAMAT ALGIPSGKHI LRRHLAEEMS TIFHPSNSTG IRESERFGTP PGRASSITRA GREDSGGNSK HRAGRLPLGK IGRGFSSKEP DFHDDYGSLE NEDCGDEDLQ GRPEQCRLEG YGSLEVTNV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Kazn
<b>Protein Names</b>	Recommended name: Kazrin
<b>Expression Region</b>	1-779
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.