



Recombinant *Propionibacterium acnes* Homoserine kinase (thrB)

Product Code	CSB-EP727996PSG-B
Abbreviation	thrB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6A8B2
Product Type	Recombinant Protein
Immunogen Species	<i>Propionibacterium acnes</i> (strain KPA171202 / DSM 16379)
Purity	≥85% (SDS-PAGE)
Sequence	MDNSVRVRVP ATSANLGPGY DCIGLALDLW DEVSVGVLDR PGVMIDVTGE GADTVPHDES HLVMATLRQG LVELGYPHPD AGLHLTAINS IPQSRGLGSS AAAVVSGLAL AWGLARPGFP LDRSALLTMA AAIEGHPDNA APAILGGAQL AWLDGEAVNH IGLTVNPSIV FRVYVPDLV PTALARQVLP EQVDRVDAVH QVLAASLLVT ALTTSPHELL AATQDWIHQY YRRALMPESA ALTDRLRGRG VATVISGAGP TVLALGSRDQ LEKVSDVDTA GFVAHDLVLG EGVHFF
Source	<i>E.coli</i>
Target Names	thrB
Protein Names	Recommended name: Homoserine kinase Short name= HK Short name= HSK EC= 2.7.1.39
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.