



Recombinant Leifsonia xyli subsp. xyli Guanylate kinase (gmk)

Product Code	CSB-EP728033LAAC-B
Abbreviation	gmk
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AF82
Product Type	Recombinant Protein
Immunogen Species	Leifsonia xyli subsp. xyli (strain CTCB07)
Purity	≥85% (SDS-PAGE)
Sequence	MEGRIRVTEF LLSVPAIGTT KTQDALRRLG IAPSKRLGGL GRHQRLRLRV FLIEREERPVR RRRKHLVLA GPTAVGKGTV SAYIRENYPE VLLSVSATTR SPRPGEVDGV NYFVDDAAF DRMIAAGDLL EHATVHNAYR YGTPRAPIEK ALDDGRSVLL EIDLQGARQV RASMPEARLI FLLPPTWEEL VRRLTGRGTE DAAEQRRLE TAKVELAAQN EFDHRVNHV VADAAREVVD LMMVRPAPVR PSTA
Source	E.coli
Target Names	gmk
Protein Names	Recommended name: Guanylate kinase EC= 2.7.4.8 Alternative name(s): GMP kinase
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.