



Recombinant Rat Serine protease 23 (Prss23)

Product Code	CSB-MP728110RA
Abbreviation	Prss23
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AY61
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	YNVPWKP TWPAYRLPIV LPQSTLKLAK PDFDAKAKLE VSSSCGPQCH KGTPLPTYEE AKQYLSYETL YANGSRTETQ VGIYILSNGE GRARSRDSEA AGKSRRKRQI YGYDGRFSIF GKDFLLNYPF STSVKLSTGC TGTLVAEKHV LTAAHCIHDG KTYVKGQTQKL RVGFLKPKYK DGAGGDNSSS SALVEKMKFQ WIRVKRTHVP KGWIKGNAND IGMDYDYALL ELKKPHKRKF MKIGVSPPAK QLPGGRIHFS GYDNDRPGNL VYRFCDVKDE TYDLLYQQCD AQPASGSGV YVRMWKRPQQ KWERKIIGIF SGHQWVDMNG SPQDFNVAVR ITPLKYAQIC YWIKNYLDG REG
Source	Mammalian cell
Target Names	Prss23
Protein Names	Recommended name: Serine protease 23 EC= 3.4.21.-
Expression Region	24-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the trypsin family of serine proteases. Mouse studies found a decrease of mRNA levels after ovulation was induced. This gene seems to be highly conserved in vertebrates and may be an important ovarian protease.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.