



Recombinant Rat THO complex subunit 6 homolog (Thoc6)

Product Code	CSB-EP728111RA-B
Abbreviation	Thoc6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AY87
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MEHATPLAVP LGQAEVFQAL QRLHMTIFSQ SVSPCGKFLA AGNNYGQIAI FSLSAALSSE AKEESKKPMV TFHAHDGPVY SMVSTDRHLL SAGDGEVKGW LWAEILKKGC KELWRRQPPY RTSLEVPEIN ALLLVPKENS LILAGGDCQL HTMDLETGTF TRALRGHTDY IHCLALRERS PEVLGGEDG AVRLWDLRIA KEVQTIEVYK HEECSRPHNG RWIGCLATDS DWMVCGGGPA LTLWHLRSST PTTVFPIRAP QKHVTFYQDL ILSAGQGCCV NHWQLSGELK AQVPGSSPGL LSLSLNQQPA APECKVLTAS GNSCRVDVFT NLGYRAFSLS F
Source	E.coli
Target Names	Thoc6
Protein Names	Recommended name: THO complex subunit 6 homolog Alternative name(s): WD repeat-containing protein 58
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.