



Recombinant Debaryomyces hansenii Signal recognition particle SEC65 subunit (SEC65)

Product Code	CSB-BP728182DIS
Abbreviation	SEC65
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BM93
Product Type	Recombinant Protein
Immunogen Species	Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)
Purity	≥85% (SDS-PAGE)
Sequence	MAKRPVLEDV LDDEIDNMDM DIAQFDPSLR TPIAPVQPGP QITRSQDQEP PLFPNFPHPE QDKSVERKDI VDPNKFSAEE RAEKKKFQII YPCYFDKSRS HKDGRRVSES KAVSNPLAKT ISDACRHYNL PVMLELDKTH PQDFGNPGRV RVLIKDNKNK GMPVDSRFKT KRAVLNIIAD YLKTHPTTLA SIGPKSGIPL PSEYETGFEP GEIPKVKGFK MNTIVPVHSE LTIKHPMTKS IYDPEPEQPQ VKAPQAPKQP KKKVMKIRG
Source	Baculovirus
Target Names	SEC65
Protein Names	Recommended name: Signal recognition particle SEC65 subunit
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.