



Recombinant *Debaryomyces hansenii* Autophagy-related protein 3 (ATG3)

Product Code	CSB-MP728203DIS
Abbreviation	ATG3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BSC4
Product Type	Recombinant Protein
Immunogen Species	<i>Debaryomyces hansenii</i> (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (<i>Torulasporea hansenii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MLRSKLSSLR EYLTPIRHTS DFTTTTGEISP EEFVEAGDYL VYKFPTWQWS SAPDKLKKDF LPPDKQYLIT KHVSSYQRAV TYLGIKSDLD EDEEELEDGW VKSHKINNDP SRLKTDASGE NSGGNTNDND DTNEINDIDE LIDENAEEQS SEDENDFEEL VETNANSNLR KYDLYITYST SYRVPKMYLV GFNSNGIPLL PKQMFEDISG DYRDKTATIE TLPVSYNTMS VSIHPCKHSS VMKVLMAHAA ASKKRENPAD LTDLAEQTGS LNLKDKVPGR DFGEDVDIEE NVPGIRVDQY LIIFLKFIAS VTPGIEYDYT MDAL
Source	Mammalian cell
Target Names	ATG3
Protein Names	Recommended name: Autophagy-related protein 3 Alternative name(s): Autophagy-related E2-like conjugation enzyme ATG3
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.