



Recombinant *Yarrowia lipolytica* Protein DML1 (DML1)

| | |
|--------------------------|---|
| Product Code | CSB-EP728260YAP-B |
| Abbreviation | DML1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6C6D9 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MREITLSFG QTAGHINADY FNSQEHYFPL ATKSTSDPTV RFRRGVATDS RTETYNPRLL TWELKGGYGA FQAFNQFYTE GDKSQPVWN EGEIKEIQE PVDKNEYQKA LDLGRENASQ LNTDTTKRWT DYNRLFHHPK TRHQLDNWLF DPDTAPQGIH RGGDQKWTGF DVGVNEWEHV LNSDKEYLDS TLRSWVEECD SLGGLNVVVD DSAWAGVAAK ILANYRDDFD AKGTVVTWSV EAKPEKKTRE TQKNAIQTTV ALSQVSSIYI PVSFPTKTP I DADFPTSPW HQAALFGLSV YEPVQLMATQ RTGGLSLDAM AALLDTTGER NIVSEVGGTL TDVDGDKSTV QYDSRLFQEL FLPSRIPPKT PHVFSDLSIA RLEDQDDDAN SYRLNQPMPE VSSQPEFTKQ FVSSKTRLAV TTKPRRTLLD MAKFVSVYER GDEREEMKQE LGDMASKYEH GWEEESDDDD DY |
| Source | E.coli |
| Target Names | DML1 |
| Protein Names | Recommended name: Protein DML1 |
| Expression Region | 1-472 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.