



Recombinant Kluyveromyces lactis Ribosome assembly protein 3 (RSA3)

Product Code	CSB-MP728316KBK
Abbreviation	RSA3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CK62
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MSEGHIAAIS KTNKRKNRRK KRRTQDVSDS SSDSSSDSEN SSVEDNDTKL ETVEKAASDV EDVTLSDIDM DKEEDEAALD DKKASRKEGS VLSSETVDKL NKVGLTTS DL TGNNGIHLGN IDLNKMSATL DESSNKLLN NKDKNGLKNQ YLSMLFESYG DDMNELRNAP DFTSKTLVML ANVLKDGGM FDLETLKTIV EDK
Source	Mammalian cell
Target Names	RSA3
Protein Names	Recommended name: Ribosome assembly protein 3
Expression Region	1-203
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.