



Recombinant *Kluyveromyces lactis* RuvB-like helicase 1 (RVB1)

Product Code	CSB-MP728346KBK
Abbreviation	RVB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CQA9
Product Type	Recombinant Protein
Immunogen Species	<i>Kluyveromyces lactis</i> (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (<i>Candida sphaerica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MVQISEVQEQ SSTAYTRTAA HTHIKGLGLD EFGVAKQVEG GFVGGQAEARE ACGVIVDLIK AKKMSGKAIL LAGGPSTGKT ALALAIQEL GPKVPCPLV GSELYSVEVK KTEALMENFR RAIGLRIKET KEVYEGEVTE LTPAEAENPL GGYGTISHV IVGLKSAKGT KTLRLDPTIY ESIQREKVSVDGVIYIESNT GAVKRVGRSD AFATEFDLEA EEYVPLPKGE VHKKKEIVQD VTLHDLVDAN ARPQGGQDVI SMMGQLMKPK KTEITEKLRH EVNKVVAKYI DQGVAELVPG VLFIDEVNML DIEIFTYLNRLAESDIAPVV VLASNRGMIT VRGTDDVVSP HGVPDLIDR LLIVRTLPLYN REEIKTIISK RAAVENLQVE DEALEFLATL GTETSLRYVL QLLSPSGIIA KIANRAEISV ADVEEAKLLF LDAKRSTKIL EQSENYL
Source	Mammalian cell
Target Names	RVB1
Protein Names	Recommended name: RuvB-like helicase 1 EC= 3.6.4.12
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.