



Recombinant Kluyveromyces lactis SWI5-dependent HO expression protein 2 (SHE2)

Product Code	CSB-EP728376KBK
Abbreviation	SHE2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CVF1
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MAYCQVVISK EILDAVLQLV STYSNYISRY VDYLNNLIAV QRRVSTLRFE RMTLIKYYVKK LRFMADVLHS WSFEGENDLI GKRLNEVIDP LGAYLIKVFE ILDLLNFYIT QPMRGETISK TLNEDLVVSE ETIVCINDSY RIYIKGIQWL VESISERNGN CPHLALIE FTRKCAIEDE VDFTASEDIL LQDVAIVEVE EEYVDLLHDW SAILVQKQTE MKELFSEDSK RWASMTKPVK K
Source	E.coli
Target Names	SHE2
Protein Names	Recommended name: SWI5-dependent HO expression protein 2
Expression Region	1-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.