



Recombinant *Erwinia carotovora* subsp. atroseptica Phosphate import ATP-binding protein PstB 2 (pstB2)

Product Code	CSB-YP728395EAAB
Abbreviation	pstB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CYN3
Product Type	Recombinant Protein
Immunogen Species	Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MSMTTETSTN KIQVRDLNFY YGKFHALKNI TLDIAANQVT AFIGPSGCGK STLLRTLNMK YQLYPEQRAE GDILLDGNNI LTDKQDIAL RAKVGMVFQK PTFPMSIYD NIAFGVRLFE KLSRADMDER VQWALTKAAL WQETKDKLHQ SGYSLSGGQQ QRLCIARGIA IRPDVLLLDE PCSALDPIS GRIEELISEL KKDYTVVIVT HNMQQAARCS DHTAFMYLGE LIEFSDTDTL FTAPRQKQTE DYITGRYG
Source	Yeast
Target Names	pstB2
Protein Names	Recommended name: Phosphate import ATP-binding protein PstB 2 EC=3.6.3.27 Alternative name(s): ABC phosphate transporter 2 Phosphate-transporting ATPase 2
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.