



Recombinant *Erwinia carotovora* subsp. atroseptica UPF0234 protein ECA1137 (ECA1137)

Product Code	CSB-MP728465EAAB
Abbreviation	ECA1137
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6D838
Product Type	Recombinant Protein
Immunogen Species	Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MPSFDIVSEI DMQEVRNAVE NATRDLSTRW DFRNVPASFE LNEKSQSIKV ASESDFQVQQ LVDILREKLV KRGIEGGALE IPEEMEHSGK TYSVEAKLKQ GIETTLAKQL IKLIKDSKLK VQAQIQGEQV RVTGKARDDL QGVMALIRGG NLGQPFQFTN FRD
Source	Mammalian cell
Target Names	ECA1137
Protein Names	Recommended name: UPF0234 protein ECA1137
Expression Region	1-163
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.