



Recombinant Arabidopsis thaliana U-box domain-containing protein 62 (PUB62)

Product Code	CSB-MP728491DOA
Abbreviation	PUB62
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DBN5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSGKGLSPAP FNAQPMIVQD ADPLRFRVGE QDPKTREFAA FIGDHRYFAA AAAAAAAAANP HPHLEFRQNF YSEKPIIGNP NDSGGSDGED DVDVEEDED DDLDGNEGDI GMNKDAGEDS VSAGAVIVVG QDNAAYYSQH FKTVEASFVS RNEESSIAAD NGCDFSGRRD LSSSSNSIE SLRTILSDPT TGSLMADAMI LPCGHTFGAG GIEQVKQMKC CCTCSQPVSE DSITPNLTLR VAVQAFCEE NSQSNHPSKR KREGFDQERR AFGVTNHSGR SRNKS NHFPF AVADRVIKIG NKRTPPRFVG REAVTTQCL NGWYVVKTLN NAESVKLQYR SLAKAPEDPS AKATPNKMVS NWL
Source	Mammalian cell
Target Names	PUB62
Protein Names	Recommended name: U-box domain-containing protein 62 EC= 6.3.2.- Alternative name(s): Plant U-box protein 62
Expression Region	1-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.