



Recombinant *Xenopus laevis* Protein FAM76B (fam76b)

Product Code	CSB-MP728511XBE
Abbreviation	fam76b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DEE7
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MASALYACTK CHQRYPFEEEL SQGQQLCKEC RIAHPIVKCT YCRSEFQQES KTNTICKKCA QNVKQFGTPK PCQYCNIAA FIGTKCQRCT NSEKKYGPPQ TCEQCKQQCA FDRKEEGRRK VDGKLLCWLC TLSYKRVLQK TKEQRKSLGS THSNSSSSSL TEKDQHHSKH HHHHHHHHRH SSSHHKISSL SPEPEQGIWK QSHKSSTVIQ NETPKKKPKL EFKPSNGDSS SINQSTDGGG TDNFVLISQL KEEVMSLKRL LQQRDQTILE KDKKLTTELKA DFQYQENNLR TKMNGMDKAH KDFVEQLQGK NRELLKQVAA LSKGKKFDKS GSILTSP
Source	Mammalian cell
Target Names	fam76b
Protein Names	Recommended name: Protein FAM76B
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.