



# Recombinant Rat Acyl-CoA-binding domain-containing protein 4 (Acbd4)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP728525RA                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Abbreviation</b>      | Acbd4                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |
| <b>Uniprot No.</b>       | Q6DGF9                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Rattus norvegicus (Rat)                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MGTEKEEPDC QKQFQAAVSV IQNLPKNGSY RPSYEEMLR FYSYYKQATAG<br>PCLVPRPGFW DFIGRYKWDA WNSLGKMSRE EAMSAYITEM KLVAQKVIDT<br>VPLGEVAEDM FG YFEPLYQV IPDMRPPET FLRRVTGWQE PAVNRDVQAA<br>PEPSHPPKEP APPSPESRLP RDL DLEV FCD SVEQLEPELV RVPVLSPVPA<br>ESELPHLHTG TGDSAQRVWA EQKEAAGREL TTRSSPESPE<br>GFGGSLMGPQ ELDRWLVGTV QAMQESMKDV HRRLQILESK<br>PQPLEQQRSP RTRPWPLGLS TPTLLFFILW PFVVQWLFRQ FRTQRR                  |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Target Names</b>      | Acbd4                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: Acyl-CoA-binding domain-containing protein 4                                                                                                                                                                                                                                                                                                                                              |
| <b>Expression Region</b> | 1-326                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Details</b>    | This gene encodes a member of the acyl-coenzyme A binding domain containing protein family. All family members contain the conserved acyl-Coenzyme A binding domain, which binds acyl-CoA thiol esters. They are thought to play roles in acyl-CoA dependent lipid metabolism. Multiple transcript variants encoding different isoforms have been found for this gene.                                      |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,                                                                                                                                                                                                                                                                                                                               |



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.