



Recombinant Human Putative histone H2B type 2-C (HIST2H2BC)

Product Code	CSB-YP728555HU
Abbreviation	HIST2H2BC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DN03
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	PEPAKFAPA PKKGSKKAVT KAQKKDGKKR KRSRKESYSI YVYKVLKRVH PDTGIWCKAM GIMNSFLNDI FERIAGEASR LAHYNKRSTI TSRRSRRPCA CCCPASWPST PCPRAPRRSP STPAPSESLP GPGARSLPPS LPPRVAGCFV SKGSFQGHLT TSVKESFLCC QSQLMFLASR LVNFRAHNT KHR
Source	Yeast
Target Names	HIST2H2BC
Protein Names	Recommended name: Putative histone H2B type 2-C Alternative name(s): Histone H2B.t Short name= H2B/t
Expression Region	2-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.