



Recombinant *Oryza sativa* subsp. *japonica* Probable protein phosphatase 2C 13 (Os02g0255100, LOC_Os02g15594)

Product Code	CSB-EP728590OFG-B
Abbreviation	Os02g0255100, LOC_Os02g15594
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6EN45
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MVCFASLRRA LPLLLRATTT TTPRFLLPRA LSGGVGGGAA VDARALLRGH SGWRGLRVAA RMMLDSSDSA AAAGQMPPQQ RAAGAVACSA QDGGAAGYAS GGWAREDGKL KCGYSSFRGK RATMEDFYDV KLTEIDGQAV SLFGVFDGHG GPRAAEYLKE NLFENLLKHP EFLTDTKLAI SETYQKTDTD FLESESNAFR DDGSTASTAV LVGGHLYVAN VGDSRAVVSK AGKAMALSED HKPNRSDERK RIENAGGVVI WAGTWRVGGV LAMSRAFGNR LLKPFVVAEP EIQEELVNEDE LECLVLASDG LWDVVENEAA VSLAKTEDLP ESVARKLTEI AYSRGSADNI TCIVVQFHHD KTE
Source	E.coli
Target Names	PP2C13
Protein Names	Recommended name: Probable protein phosphatase 2C 13 Short name= OsPP2C13 EC= 3.1.3.16
Expression Region	1-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.